

KLKKTETQ

Chemical Properties

CAS No. : 1054611-00-5


Formula: C42H78N12O14

Molecular Weight: 975.14

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	KLKKTETQ is a bioactive peptide recognized for its anti-aging properties, commonly utilized as an ingredient in cosmetics [1].
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0255 mL	5.1275 mL	10.2549 mL
5 mM	0.2051 mL	1.0255 mL	2.051 mL
10 mM	0.1025 mL	0.5127 mL	1.0255 mL
50 mM	0.0205 mL	0.1025 mL	0.2051 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481