

JAK-IN-33

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	JAK-IN-33 (Compound 3 (R)) is a selective inhibitor of the Janus kinase (JAK) family [1].
Targets(IC50)	JAK
In vitro	JAK-IN-33 demonstrates effective utilization in the bloodstream [1].

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481