


Isohericerin

Chemical Properties

CAS No. : 140381-53-9
 Formula: C₂₇H₃₃NO₃
 Molecular Weight: 419.56
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Isohericerin is a natural isoindolinone-derived product.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3834 mL	11.9172 mL	23.8345 mL
5 mM	0.4767 mL	2.3834 mL	4.7669 mL
10 mM	0.2383 mL	1.1917 mL	2.3834 mL
50 mM	0.0477 mL	0.2383 mL	0.4767 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481