

IRGRILHDGAYSLTLQGLGIH TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	IRGRILHDGAYSLTLQGLGIH TFA, a multi-epitope class II rat HER2/neu peptide, is utilized in the study of class II HER2-DC1 vaccines [1].
-------------	---

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481