

IRBP derived peptide, R16 acetate

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

IRBP-derived peptide R16, an active biological peptide, is a segment of the interphotoreceptor retinoid-binding protein (IRBP) known to induce experimental autoimmune uveitis (EAU) in susceptible animal strains.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481