


Interiotherin D

Chemical Properties

CAS No. : 460090-66-8
 Formula: C₂₆H₂₆O₈
 Molecular Weight: 466.48
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Interiotherin D, a lignan, can be isolated from Kadsura interior [1].
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1437 mL	10.7186 mL	21.4371 mL
5 mM	0.4287 mL	2.1437 mL	4.2874 mL
10 mM	0.2144 mL	1.0719 mL	2.1437 mL
50 mM	0.0429 mL	0.2144 mL	0.4287 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481