

## Imbotolimod

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Imbotolimod, a humanized monoclonal antibody of the immunoglobulin G1-kappa class, exhibits both anti-ERBB2 and antineoplastic activities. It is a derivative of telratolimod [1].
Targets(IC50)	EGFR,Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481