

Idarucizumab

Chemical Properties

CAS No. : 1362509-93-0

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Idarucizumab, a humanized monoclonal antibody fragment, serves as the first reversal agent for direct oral anticoagulants (DOACs). It specifically neutralizes the effects of oral direct thrombin inhibitors, thereby restoring hemostasis [1].
Targets(IC50)	Thrombin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481