


ICMT-IN-6

## Chemical Properties

CAS No. : 1313602-75-3  
 Formula: C<sub>23</sub>H<sub>31</sub>NO<sub>2</sub>  
 Molecular Weight: 353.5  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	ICMT-IN-6 (compound 29) is an inhibitor of isoprenylcysteine carboxyl methyltransferase (ICMT) with an IC <sub>50</sub> value of 0.09 μM [1].
Targets(IC <sub>50</sub> )	Others,Transferase

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8289 mL	14.1443 mL	28.2885 mL
5 mM	0.5658 mL	2.8289 mL	5.6577 mL
10 mM	0.2829 mL	1.4144 mL	2.8289 mL
50 mM	0.0566 mL	0.2829 mL	0.5658 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481