

ICMT-IN-47

Chemical Properties


CAS No. : 1313602-72-0

Formula: C₂₅H₃₅NO

Molecular Weight: 365.55

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	ICMT-IN-47 (compound 26) is an ICMT inhibitor with an IC ₅₀ value of 0.76 μM [1].
Targets(IC ₅₀)	Others,Transferase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7356 mL	13.678 mL	27.356 mL
5 mM	0.5471 mL	2.7356 mL	5.4712 mL
10 mM	0.2736 mL	1.3678 mL	2.7356 mL
50 mM	0.0547 mL	0.2736 mL	0.5471 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481