

Human IgG4 (S228P) lambda1, Isotype Control

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Human IgG4 (S228P) lambda1 Isotype Control is a chimeric monoclonal antibody of human-mouse origin, serving as an isotype control specifically for the human IgG4 (S228P) λ1 antibody.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481