

## Human IgG1 lambda1, Isotype Control

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Human IgG1 lambda1, Isotype Control is a human IgG1 subtype monoclonal antibody that can be used as a negative control or reference antibody to assess the specificity of primary antibodies in ELISA, flow cytometry, and immunohistochemistry.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481