

## Human IgG1 (D270N) kappa, Isotype Control

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

-20°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Human IgG1 (D270N) kappa Isotype Control is a chimeric human-mouse monoclonal antibody that serves as an isotype control for the human IgG1 (D270N)κ antibody.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481