

HSP70/SIRT2-IN-1

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	HSP70/SIRT2-IN-1 (Compound 2a) serves as a dual inhibitor targeting both SIRT2 and HSP70, exhibiting an IC50 of 17.3±2.0 µM against SIRT2. This compound possesses antitumor activity [1].
Targets(IC50)	HSP,Sirtuin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481