


Hosenkoside N

Chemical Properties

CAS No. :	156765-13-8
Formula:	C42H72O15
Molecular Weight:	817.01
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Hosenkoside N, a baccharane glycoside, has been isolated from Impatiens balsamina [1].
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.224 mL	6.1199 mL	12.2398 mL
5 mM	0.2448 mL	1.224 mL	2.448 mL
10 mM	0.1224 mL	0.612 mL	1.224 mL
50 mM	0.0245 mL	0.1224 mL	0.2448 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481