

HIF-1 alpha (556-574) (TFA)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	HIF-1 alpha (556-574) TFA is a 19-residue fragment of hypoxia-inducible factor-1 (HIF-1), the master regulator of oxygen homeostasis [1].
Targets(IC50)	HIF/HIF Prolyl-Hydroxylase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481