

Hexapeptide-42

Chemical Properties

CAS No. : 197771-66-7


Formula: C₂₅H₄₆N₁₀O₈S

Molecular Weight: 646.76

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

Hexapeptide-42 is a bioactive peptide employed as a cosmetic ingredient [1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5462 mL	7.7308 mL	15.4617 mL
5 mM	0.3092 mL	1.5462 mL	3.0923 mL
10 mM	0.1546 mL	0.7731 mL	1.5462 mL
50 mM	0.0309 mL	0.1546 mL	0.3092 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481