


## Hexapeptide-3

## Chemical Properties

CAS No. :	1205679-02-2
Formula:	C32H58N14O11S
Molecular Weight:	846.96
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Hexapeptide-3 is a bioactive peptide known for its skin-repairing effects and is documented as a cosmetic ingredient [1].
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1807 mL	5.9035 mL	11.8069 mL
5 mM	0.2361 mL	1.1807 mL	2.3614 mL
10 mM	0.1181 mL	0.5903 mL	1.1807 mL
50 mM	0.0236 mL	0.1181 mL	0.2361 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481