

## H-Asp(AMC)-OH

## Chemical Properties

CAS No. : 133628-73-6


Formula: C<sub>14</sub>H<sub>14</sub>N<sub>2</sub>O<sub>5</sub>

Molecular Weight: 290.27

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	H-Asp(AMC)-OH, an amino acid derivative and fluorescent dye, does not inhibit glycine transport at a concentration of 0.25 mM [1].
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4451 mL	17.2253 mL	34.4507 mL
5 mM	0.689 mL	3.4451 mL	6.8901 mL
10 mM	0.3445 mL	1.7225 mL	3.4451 mL
50 mM	0.0689 mL	0.3445 mL	0.689 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481