

Gotistobart

Chemical Properties

CAS No. : 2226344-78-9

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Gotistobart (ONC-392; BNT 316) is a humanized monoclonal antibody targeting CTLA-4/CD152 that selectively depletes regulatory T cells (Tregs) in the tumor microenvironment and enhances immunotherapy efficacy, applicable for studying non-small cell lung cancer (NSCLC).
Targets(IC50)	Immunology/Inflammation related

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481