

Golodacimab

Chemical Properties

CAS No. : 2418540-63-1

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Golodacimab (MEDI6570) is a fully human monoclonal antibody that specifically targets the lectin-like oxidized low-density lipoprotein receptor-1 (LOX-1), with research indicating potential therapeutic applications for mitigating complications associated with type 2 diabetes mellitus and for the prophylaxis of hereditary angioedema (HAE) attacks.
Targets(IC50)	Lipoxygenase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481