


Glyurallin B

Chemical Properties

CAS No. :	199331-53-8
Formula:	C ₂₅ H ₂₆ O ₆
Molecular Weight:	422.47
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Glyurallin B, a flavonoid derived from Glycyrrhiza inflata and Glycyrrhiza uralensis (licorice), exhibits ABTS+ radical scavenging activity and lipid peroxidation inhibition, with EC ₅₀ values of 11.9 ± 0.58 μM and 15.3 ± 1.26 μM, respectively [1].
Targets(IC ₅₀)	Reactive Oxygen Species
In vitro	Glyurallin B, at a concentration of 10 μM, exhibits significant cytotoxicity against RAW 264.7 cells, with cell viability measured at 26.3 ± 5.3% [1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.367 mL	11.8352 mL	23.6703 mL
5 mM	0.4734 mL	2.367 mL	4.7341 mL
10 mM	0.2367 mL	1.1835 mL	2.367 mL
50 mM	0.0473 mL	0.2367 mL	0.4734 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481