

GluN1(359-378)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Keep away from moisture

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	GluN1 (359-378), an antibody targeting the N-methyl-D-aspartate receptor (NMDAR), can cross the blood-brain barrier and is used in research on immune-based therapies for anti-NMDAR encephalitis [1].
Targets(IC50)	iGluR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481