


## Ganoderic acid Me

## Chemical Properties

CAS No. : 108026-93-3  
 Formula: C<sub>34</sub>H<sub>50</sub>O<sub>6</sub>  
 Molecular Weight: 554.76  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Ganoderic acid Me, a triterpenoid present in Ganoderma lucidum [1], exhibits various biological activities.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8026 mL	9.0129 mL	18.0258 mL
5 mM	0.3605 mL	1.8026 mL	3.6052 mL
10 mM	0.1803 mL	0.9013 mL	1.8026 mL
50 mM	0.0361 mL	0.1803 mL	0.3605 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481