


Galanin-B2

Chemical Properties

CAS No. :	944087-28-9
Formula:	C104H174N24O22
Molecular Weight:	2112.64
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Galanin-B2 (NAX-5055) demonstrates activity in epilepsy models and possesses the ability to penetrate the blood-brain barrier [1].
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.4733 mL	2.3667 mL	4.7334 mL
5 mM	0.0947 mL	0.4733 mL	0.9467 mL
10 mM	0.0473 mL	0.2367 mL	0.4733 mL
50 mM	0.0095 mL	0.0473 mL	0.0947 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481