


## G-{d-Arg}-GDSPASSK

## Chemical Properties

CAS No. :	133640-14-9
Formula:	C37H64N14O16
Molecular Weight:	960.99
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	G-{d-Arg}-GDSPASSK is a polypeptide that inhibits cell adhesion, facilitating wound healing and promoting tissue repair [1].
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0406 mL	5.203 mL	10.4059 mL
5 mM	0.2081 mL	1.0406 mL	2.0812 mL
10 mM	0.1041 mL	0.5203 mL	1.0406 mL
50 mM	0.0208 mL	0.1041 mL	0.2081 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481