

## Firefly Luciferase circRNA

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Keep away from direct sunlight

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Firefly Luciferase circRNA, originally isolated from Photinus pyralis (Firefly), expresses a luciferase protein and serves as an optimal marker for monitoring and optimizing transfection efficiency; it is also recommended as a positive control for circular RNA transfections.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481