

FAK-IN-12

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	FAK-IN-12 (Compound 125) is a selective FAK inhibitor with an IC50 of 47 nM, effectively suppressing proliferation of MGC-803, HCT-116, and KYSE30 cells with IC50 values of 0.24, 0.45, and 0.44 $\mu$ M, respectively. Additionally, FAK-IN-12 promotes apoptosis and cellular senescence [1].
Targets(IC50)	FAK

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481