

Etigilimab

Chemical Properties

CAS No. : 2044984-83-8

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Etigilimab represents a highly specialized humanized IgG1 monoclonal antibody specifically engineered to target the T-cell immunoreceptor with Ig and ITIM domains (TIGIT) with exceptional binding affinity and selectivity. It effectively facilitates the measurable blockade of the TIGIT/CD155 (PVR) inhibitory axis across various preclinical experimental models to evaluate the restoration of T-cell and NK-cell effector functions during strictly monitored laboratory observation periods and functional immune activation assays to ensure high experimental sensitivity.
Targets(IC50)	Immunology/Inflammation related

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481