

Epoetin beta

Chemical Properties

CAS No. : 122312-54-3

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Epoetin beta (rhEPO), a recombinant erythropoietin, is crucial for sustaining erythropoiesis and is utilized in anemia research [1].
Targets(IC50)	Others
In vivo	Epoetin beta (rhEPO) administered subcutaneously at a dosage of 1000 IU/kg three times weekly for 14 days prevents anemia and increases radiosensitivity in rapidly proliferating DS-sarcomas [2].

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481