

Endotrophin (Mus musculus)

Chemical Properties

CAS No. : 1678414-54-4

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Endotrophin (Mus musculus), a cleavage product of collagen 6 (Col6), induces upregulation of profibrotic and proinflammatory genes [1].

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481