

## Eglatoprutug

### Chemical Properties

CAS No. : 2518200-21-8

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Eglatoprutug is a humanized immunoglobulin G1-kappa monoclonal antibody targeting the c-Kit receptor [1].
Targets(IC50)	c-Kit

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481