

Efaprimin alfa

Chemical Properties

CAS No. : 2102200-64-4

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | Efaprimin alfa (OMP-336B11) is a human monoclonal antibody targeting TNFRSF18, functioning as a GITR ligand-Fc fusion protein [1]. |
| Targets(IC50) | TNF |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481