

## Efanesoctocog alfa

## Chemical Properties

CAS No. : 2252477-42-0

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Efanesoctocog alfa (BIVV001) is a humanized antibody with continuous procoagulant cofactor activity (CPCA) [1].
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481