

## Edobacomab

## Chemical Properties

CAS No. : 141410-98-2  
Formula:  
Molecular Weight:  
Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Edobacomab (E5), a monoclonal IgM antibody targeting endotoxin, is utilized in the study of sepsis stemming from Gram-negative bacterial infections [1].
Targets(IC50)	Others
In vitro	Prolgolimab has the capacity to bind to endotoxin, neutralizing it, which consequently attenuates the inflammatory response and leads to improved patient outcomes [2].

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481