

Ecleralimab

Chemical Properties

CAS No. : 2415207-91-7

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Ecleralimab is a potent inhaled antibody fragment (scFv) targeting TSLP, acting by blocking the binding of TSLP to its receptor complex. The compound is an essential tool for investigating the role of airway epithelial-derived cytokines in asthma pathogenesis, aiming to alleviate Type 2 inflammatory responses by inhibiting the upstream immune cascade.
Targets(IC50)	Interleukin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481