


DGK ζ -IN-3

Chemical Properties

CAS No. :	2732901-70-9
Formula:	C ₂₀ H ₁₆ ClF ₃ N ₄ O ₃ S
Molecular Weight:	484.88
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	DGK ζ -IN-3 is an inhibitor of diacylglycerol kinase zeta (DGK- ζ), employed in the treatment of both liquid and solid cancers, as well as other diseases regulated by DGK- ζ [1].
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0624 mL	10.3118 mL	20.6237 mL
5 mM	0.4125 mL	2.0624 mL	4.1247 mL
10 mM	0.2062 mL	1.0312 mL	2.0624 mL
50 mM	0.0412 mL	0.2062 mL	0.4125 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481