

Delpacibart

Chemical Properties

CAS No. : 2653183-66-3

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Delpacibart is a humanized monoclonal antibody targeting TFRC (Transferrin Receptor), applicable for Duchenne muscular dystrophy (DM) research and commonly used in synthesizing AOCs (Antibody Oligonucleotide Conjugates) such as Delpacibarat etedesiran.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481