

DEALA-Hyp-YIPD

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

DEALA-Hyp-YIPD, an HIF-1 α peptide inhibitor, impedes the VHL/HIF-1 α interaction with an half-maximal inhibitory concentration (IC₅₀) of 0.91 μ M and a dissociation constant (K_d) of 180 nM [1].

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481