

DDa-1


## Chemical Properties

CAS No. : 3036944-98-3

Formula: C60H77Cl2N13O3S

Molecular Weight: 1131.31

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.



## Biological Description

Description	DDa-1 is a potent (kinase degrader) [1].
Targets(IC50)	Others,BTK,PROTACs

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8839 mL	4.4197 mL	8.8393 mL
5 mM	0.1768 mL	0.8839 mL	1.7679 mL
10 mM	0.0884 mL	0.442 mL	0.8839 mL
50 mM	0.0177 mL	0.0884 mL	0.1768 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481