

DCAF1 binder 2

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	DCAF1 Binder 2, an E3 ligase ligand, is implicated in BRD9, kinase, and selective Bruton's tyrosine kinase (BTK) degradation.
Targets(IC50)	Others,Epigenetic Reader Domain,BTK,Ligands for E3 Ligase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481