

## Dazodalibep

## Chemical Properties

CAS No. : 2245953-10-8

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Dazodalibep (MEDI 4920; VIB 4920) is a monoclonal antibody that specifically targets CD40LG/TNFSF5 and is fused to human serum albumin (ALB/HSA) [1].
Targets(IC50)	TNF

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481