

Davutamig

Chemical Properties

CAS No. : 2648058-48-2

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Davutamig (REGN-5093) is a humanized IgG4-kappa anti-MET monoclonal antibody that targets two distinct, nonoverlapping epitopes on the MET receptor and functions as an antineoplastic agent [1].
Targets(IC50)	c-Met/HGFR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481