

## Dapirolizumab pegol

## Chemical Properties

CAS No. : 1416147-64-2

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Dapirolizumab pegol (CDP 7657), a monoclonal antibody of the Fab'-IgG-κ class, targets CD40LG/TNFSF5 and is conjugated with a pegol moiety. The antibody was initially cloned from the Norway rat ( <i>Rattus norvegicus</i> ).
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481