

Cy3-PEG3-TCO4

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy3-PEG3-TCO4, a dye derivative with three polyethylene glycol (PEG) units linked to Cyanine 3 (Cy3), possesses a trans-cyclooctene (TCO) group capable of undergoing the inverse electron demand Diels-Alder reaction (IEDDA) with tetrazine-functionalized molecules.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481