

## Cy3-PEG3-TCO

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Cy3-PEG3-TCO, a Cyanine 3 (Cy3) dye derivative featuring three polyethylene glycol (PEG) units, utilizes its TCO moiety to engage in inverse electron demand Diels-Alder reactions (iEDDA) with tetrazine-functionalized molecules.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481