


## Condurango glycoside E3

## Chemical Properties

CAS No. :	115784-10-6
Formula:	C66H98O26
Molecular Weight:	1307.47
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Condurango Glycoside E3 (Compound 5), a pregnane glycoside, is isolated from the dried bark of Marsdenia cundurango [1].
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7648 mL	3.8242 mL	7.6484 mL
5 mM	0.153 mL	0.7648 mL	1.5297 mL
10 mM	0.0765 mL	0.3824 mL	0.7648 mL
50 mM	0.0153 mL	0.0765 mL	0.153 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481