


Codactide

Chemical Properties

CAS No. :	22572-04-9
Formula:	C101H158N30O23S
Molecular Weight:	2192.59
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Codactide is a synthetic octadecapeptide exhibiting corticotropin-like activity [1].
Targets(IC50)	Melanocortin Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.4561 mL	2.2804 mL	4.5608 mL
5 mM	0.0912 mL	0.4561 mL	0.9122 mL
10 mM	0.0456 mL	0.228 mL	0.4561 mL
50 mM	0.0091 mL	0.0456 mL	0.0912 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481