


Cinnamtannin D2

Chemical Properties

CAS No. :	97233-47-1
Formula:	C60H48O24
Molecular Weight:	1153.01
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Cinnamtannin D2 (compound 13), a tannin derived from <i>Vaccinium vitis-idaea</i> L., boasts isolation from the aforementioned source [1].
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8673 mL	4.3365 mL	8.673 mL
5 mM	0.1735 mL	0.8673 mL	1.7346 mL
10 mM	0.0867 mL	0.4336 mL	0.8673 mL
50 mM	0.0173 mL	0.0867 mL	0.1735 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481