


Caffeoyl-CoA

Chemical Properties

CAS No. :	53034-79-0
Formula:	C30H42N7O19P3S
Molecular Weight:	929.68
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Caffeoyl-CoA can be isolated from the leaf tissues of Gramineae, exhibiting a Km value of 1.45 μ M [1].
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0756 mL	5.3782 mL	10.7564 mL
5 mM	0.2151 mL	1.0756 mL	2.1513 mL
10 mM	0.1076 mL	0.5378 mL	1.0756 mL
50 mM	0.0215 mL	0.1076 mL	0.2151 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481